



### Shandon Hude Dental Material Corporation

Add / No. 68 Shanhai Road, Donggang District, Rizhao City, Shandong Province,  
276800, P.R. China.  
Tel / +86 (633) 2277268    marketing@hugedental.com    www.hugedental.com

- Facebook Huge Dental
- Instagram Huge Dental
- Youtube Huge Dental

v1.1-202506

**HUGE** Your Partner  
in Prosthetic  
Dentistry

# CLINICAL SOLUTION

Redefine the fully-modernized dental clinics.





# Company Profile

HUGE Dental that integrates research, development, manufacture and sales of dental materials, was founded in 2006 in China. There are three R&D centers functioning independently in Rizhao, Shanghai and Suzhou.

Implementing strict QC standards for dental industry, HUGE Dental is one of top dental material manufacturers in China, and has passed the MDR CE(EU), FDA(USA), MDL(Canada), ISO13485 and other quality system certifications and product registrations. Its dental consumables have been exported to more than 100 countries/ regions all over the world.

HUGE Dental will continue to take independent innovation as the core motivity, to provide customers with integrated dental solutions.

 **EU-MDR FDA**  
Certified

**33,000+m<sup>2</sup>**  
Manufacturing and  
research plant

**700+**  
Employees

**100+**  
Sale Countries

**50,000,000+**  
Product Users





# Development History

## 2006-2009

Established Shandong Huge Dental Material Corporation.  
HUGE CAD/CAM tooth molding center was established, 8 series of teeth with different mould design that can be applied to all races worldwide were launched.

## 2010-2016

Achieved mass production of elastomeric impression materials.  
Began developing clinical dental products.

## 2017-2019

Our clear aligner system received Class II medical device certification and entered commercial production.  
Our Light Cure Dental Adhesive and Dual Cure Resin Cement were certified as Class III medical devices.

## 2020-2021

Both our Light Cure Dental Adhesive and Dual Cure Resin Cement received CE certification and FDA approval.  
Our Light Cure Veneer Cement obtained Class III medical device certification  
Our Universal Composite Restorative was approved by the FDA.  
HUGE DENTAL USA LLC was founded in California, USA.

## 2022

Joined the Innovation Center of the “National Engineering Research Center of Oral Biomaterials and Digital Medical Devices” .  
Our Glass Ionomer Cement was certified in Switzerland.  
Our Flowable Composite Restorative received FDA approval.

## 2023

Twelve of our products received initial registration certificates from the UK MHRA.  
Shandong Huge Dental Material Corporation was recognized as one of the “Top 100 Most Innovative Private Enterprises in Shandong Province” .  
Our Universal Composite Restorative received Class III certification from the NMPA.  
Our Glass Ionomer Cement, Light Cure Dental Adhesive, and Dual Cure Resin Cement were registered with the UK MHRA.  
Our Enamel Coating Resin was approved by the FDA.

## 2024 till Now


Established PT Huge Dental Indonesia to further expand our presence in Southeast Asia.  
Since 2024, HUGE Dental has launched a comprehensive range of digital dental equipment, offering a one-stop solution for digital dental labs and chairside applications.

# HUGE E-shop is Live!

E-Shop Europe:

[store.hugedental.com](https://store.hugedental.com)


- Multi language
- Wide range of dental products
- Exclusive Offers



E-Shop USA & Canada:

[store.hugedentalusa.com](https://store.hugedentalusa.com)

- Special Offer
- Loyalty Points and Rewards
- Fast Shipping








# HUGE All-In-One Clinical Solutions

As an innovator and reliable clinical solution provider, HUGE is committed to offering fully aligned clinical product systems work together leading to complete solutions, speeding up your path to practice success.

HUGE provides comprehensive clinical solutions, including TopCEM indirect restorative solutions, TrusFIL direct aesthetic restorative solutions, TopCEM veneer aesthetic solutions, implant restorative solutions, and VinciSmile orthodontic solutions. State-of-the-art materials, tools, and techniques which are all included in clinical solutions, result in streamlined procedures, enhanced practice efficiency and better treatment results, which in turn increase patient satisfaction and bring you to dental success.







## > Invisible Orthodontic Solution

· Digital Invisible Orthodontic Solution	03
--	----

## > Dental Restorative Materials

### 01 Pre-treatment

· HugeBond Universal FliPro Light Cure Dental Adhesive	09
· P-Etchant Phosphoric Acid Etching Gel	11
· HF-Etchant Hydrofluoric Acid Etching Gel	12
· TopCEM Ceramic Primer Ceramic Coupling Agent	13

### 02 Dental Cements

· TopCEM Dual Cure Resin Cement	15
· TopCEM Elite System-Universal Adhesive Cementation Solution	17
· TopCEM-Veneer Light Cure Veneer Cement	19
· TopCEM Veneer System-Veneer Cementation Solution	21
· TopCEM RMGI Resin Modified Glass Ionomer Cement-Luting	23
· TopCEM Vigor SA Self Adhesive Resin Cement	25

### 03 Dental Composite Restorative

· TrusFIL Universal Composite Restorative	27
· TrusFIL-Flow Flowable Composite Restorative	29
· TrusFIL X-Blend Universal Composite Restorative	31
· TrusFIL Restoration System-Universal Composite Restoration Solution	33

### 04 Core Build-up

· UltraCore Core Build-up Material	35
------------------------------------	----

### 05 Cavity Prevention

· Proseal F Enamel Coating Resin	37
· Pro Shield Flour Protector	39
· Pediatric Dental Solution	41

## > Dental Impression Materials

· PERFIT Putty	45
· PERFIT Heavy Body	47
· PERFIT Regular Body	49
· PERFIT Light Body	51
· PERFIT Super Light Body	53
· Impression Kit	57
· PERFIT A-Silicone for Bite Registration	59
· PERFIT Tray Adhesive	61

## > Glass Ionomer Cement

· Glass Ionomer Cement Luting I	75
· Glass Ionomer Cement Filling I	77
· Glass Ionomer Cement LuFill HS	79
· Glass Ionomer Cement Capsule	81

## > Intraoral Scanner

· Intraoral Scanner	85
---------------------	----

Your partner in  
 prosthetic dentistry



# ORTHODONTIC APPLIANCE





# Why VinciSmile Clear Aligner?

- 01

**Innovative Double-Aligner System**  
To ensure the ideal treatment outcomes, VinciSmile takes the soft aligner to start gently and uses the thicker aligner to control dental root movement precisely.
- 02

**Wide Range of Treatable Cases**  
Suitable to correct various malocclusions from mild to complicated cases.
- 03

**Expert Orthodontist-Led Team**  
Designed by the expert VinciSmile orthodontist-led team with scientific and accurate tooth movement, the 3D plan satisfies the requests from doctors and patients.
- 04


**Intelligent Software & System**  
Self-built system and portal providing the ortho system & dealer system. Convenient multi-interface 3D viewers, including Web, Studio, and APP. Intelligent 3D treatment plan sharing-link function, conducive to effective communication.

- 05

**Short Turnaround Time**  
VinciSmile pursues the high-efficient 1+4 mode, 1 day for 3D plan design, and 4 days for aligner fabrication.
- 06

**Professional Training Supports**  
Online and offline training courses offered to doctors for exploring and developing their potential in better medical practice.

Your Smile.  
Our Mission.



## VinciSmile® Classic

Supported by HUGE professional orthodontist team, HUGE Classic aligners are made of the high-quality biomedical polymer material to give the new smiles.



## VinciSmile® Premium

With a soft aligner and a hard aligner for each step, HUGE Premium helps speed up complex tooth movement, even more comfortable.



## VinciRetainer

Make the teeth retain their proper alignment, keeping the teeth at the ideal position, minimising any relapse.







## Self-developed Intelligent Software

Leveraging research from the National University of Singapore and the Singapore Bioimaging Center, our software benefits from highlighted features to enhance treatment precision.



### Multi Language

Available in English, Spanish, and Japanese.



### Ceph System

Incorporates 57 landmarks, 125 measurement items, and 14 analysis method.



### Quick Demo

Quickly generates simulated results of intra-oral models after treatment.



### Case Library

Covers a wide variety of successful cases for doctors to study.

## Orthodontist-Led Team

Ensure more accurate treatment plans.

Accurate, reliable and effective VinciSmile® 3D treatment plan.


2 Orthodontists with master degree



25+ 3D scheme designers



40+ Professional orthodontic team



## Process of Clinical Application

3 days for 3D treatment plan design and 4 days for appliance manufacturing.







# RESTORATIVE MATERIAL





# HugeBond Universal

## FliPro

Light Cure Dental Adhesive

A leading innovation meeting all your expectations for a universal bonding.



“EASY TO USE WITH MANY INDICATIONS.”

Indication

- Direct bonding under a light-cured composite or compomer fillings.
- Repair of dental restorations of ceramic, composites and metal alloys.
- Indirect dental restorations cementation with resin cement.
- Seal hypersensitive areas of teeth.

8<sup>TH</sup> GENERATION

Universal Application

- Contains **MDP monomer** which enhances chemical adhesion to direct and indirect restorations.
- Suitable for various substrates.
- Compatible with all etching techniques, with cement of light-cured, self-cured and dual-cured.

Reliable Bonding

- Strong bonding strength no lower than 20MPa.
- Low film thickness around 3µm.
- Ultra-mild pH(PH=2.7)provides excellentdentin bond strengths and enamel bonds.


Excellent Operative Experience

- Low occurence of microleakage prevents secondary caries.
- Flip-open design enables easy application with only one hand.
- Suction technique delivers drop control result, economizes on the output of HugeBond.


Evaluators' Comments

- “The ability to use a single bonding material for all styles and types of composites is very nice.”
- “The viscosity was thicker than most agents. I felt I had great coverage.”
- “The bottle design made it easy to get the amount you needed. Some bottle designs dispense more than what you need or want.”
- “Great consistency and clinical usability.”

DENTAL ADVISOR  
Biomaterials Test Results:  
Shear Bond Strength  
Self-Etch Mode



Evaluation Summary:





# P-Etchant

## Phosphoric Acid Etching Gel

37% Phosphoric acid etchant with ideal consistency for accurate placement

### Indication


- Used for etching of tooth enamel and/or dentin.
- Used for cleaning of restoration surfaces prior to repairing or cementation.

### Advantages

- Eliminate the surface dirt and fouling from enamel and dentin.
- Roughen the enamel and dentin surface effectively.
- Fully open dentine tubules, which is easier for adhesive to penetrate and form resin protrusion.
- Bright blue color can be easily differentiated from the tooth.
- Consistent material thickness ensures simple and effective application.
- It rinses perfectly, as it rinses off easily without leaving residue.
- Ideal viscosity, which makes precision placement easier and superior control.

### Packaging

3ml\*1 Syringe  
10 Dispensing Tips



# HF-Etchant

## Hydrofluoric Acid Etching Gel

Reliable hydrofluoric acid etchant of 9% concentration

### Indication

- Used for etching various types of resin bonded porcelain/ceramic restorations before the bonding and cementation procedure.

### Advantages

- Used for etching pretreatment of restoration adhesion surface made of glass ceramics.
- It makes the restoration surface more rough and improved adhesion.
- It enhances mechanical interlocking force.

### Packaging

1.2ml\*1 Syringe  
10 Dispensing Tips





# TopCEM Ceramic Primer

## Ceramic Coupling Agent

A Universal primer for conditioning all common restorative materials.

### Indication


- Ceramic Coupling Agent serves as priming agent and is used to create a durable adhesion between luting composites and glass/oxide ceramic, metal, composite and fiber-reinforced composite restorations.

### Advantages

- Quick solvent evaporation minimizes waiting time.
- Serving as a universal primer, suitable for multiple restoration interfaces.
- Simple operation, easy to apply, compatible with light cured or dual cured adhesive resin cements.
- Significantly improve the bonding strength between the resin cement and the restoration.

### Packaging

5ml\*1 Bottle



\* Data from HUGE Dental lab.

Around  
**25MPa**

High Bonding Strength to  
Glass Ceramics and Zirconia





# TopCEM

## Dual Cure Resin Cement

An all-rounder for all types of restorations and materials cementation.


### Indication

- Cementation of all sorts of dental restorations including crowns, bridges, inlays/onlays, veneers, dental posts and other restorations made from materials of metals/alloys, metal-ceramic, all-ceramic and/or porcelain, composites, and their combinations.

### Advantages

- Dual cure offers flexible options.
- Sufficient working time and fast setting time.
- High flexural strength over 120MPa delivers a long-lasting adhesion.
- Low film thickness about 13µm.
- Low solubility and water absorption.
- Impressive radiopacity over 400% enables easier diagnosis, treatment, and monitoring of dental restorations at recall.
- Innovative PTHP technology mitigates the high temperature sensitivity, addressing potential premature curing issues during transportation.
- PTHP redox technology remarkably enhances the color stability, reduces color shifts before and after curing.


### Shade



### Technical Parameter

Flexural Strength/ MPa	123MPa
Film thickness/µm	13µm
Water Absorption/ µg/mm³	16µg/mm³
Solubility/ µg/mm³	0.6µg/mm³
Working Time	≥2min30sec
Curing Time	2min40sec~ 3min10sec

### Bonding Strength(MPa)



### Packaging

8g\*1 Syringe





# TopCEM Elite System

## Universal Adhesive Cementation Solution

### Indication

TopCEM Elite System is designed for cementation of all sorts of dental restorations including crowns, bridges, inlays/onlays, veneers, dental posts and other restorations, which are made of metals/alloy, glass ceramic, oxide ceramic (zirconia), metal-ceramic, all-ceramic, composites, and their combinations.

### Packaging


6ml*1 bottle	HugeBond Universal FliPro Light Cure Dental Adhesive
8g*1 syringe	TopCEM Dual Cure Resin Cement
5ml*1 bottle	TopCEM Ceramic Primer Ceramic Coupling Agent
3ml*1 syringe	P-Etchant Phosphoric Acid Etching Gel
Accessories:	Intraoral mixers, dispensing tips, mixing pad, micro-brushes

### Restoration Pre-treatment

• Porcelain/ glass ceramics restorations




Apply HF-Etchant to the restoration surface. Wait for 60sec.



Rinse thoroughly with water for 10sec.



Dry with oil-free air.



Apply TopCEM Ceramic Primer to restoration inner surface. Wait for 60sec.



Dry the primed surface. \*DO NOT TOUCH THE SURFACE WITH YOUR HANDS.

• Zirconia/ metal restorations



Use aluminum oxide to sandblast bonding surfaces of restorations.



Rinse with water and dry with oil-free air.



Apply HugeBond Universal FliPro and allow it to react for 20 sec.



Dry the restorations with air.

• Composite resin restorations

Follow the same procedure as "Zirconia/ metal restorations".

### Tooth Pre-treatment



Clean the tooth to be cemented, rinse with water.



Gently dry the tooth surface with oil-free air or cotton pellets. Be careful not to over dry.

[OPTIONAL]

**A**

Enamel: 20-30s  
Dentin: Max 15s

**B**

10s

**C**


Gently dry the tooth surface with air.

### Bonding Agent Applying

Note: Do not light cure bonding agent on the indirect restoration.



Dry the tooth surface with oil-free air for 2-3sec. Be careful not to over dry.



Apply HugeBond Universal FliPro evenly in two coats for 20 sec.




Dry the tooth surface with oil-free air for 10-15sec to remove the solvent.



Light cure the tooth surface for 10sec.

### Restoration Cementation

• Crowns, bridges, inlays, onlays, veneers



Place the 1:1 mixed TopCEM onto the primed surface of the restoration.



Seat the restoration on the tooth in position.



Light cure the excess cement in the margin for 2sec, remove and clean up.



Light cure each side of the restoration for 10sec, or wait for 4.5 min for self curing.



Finish and polish, adjust the occlusal.





## TopCEM Veneer Cement

### Light Cure Veneer Cement

Your best option which brings to optimal outcomes for prosthetic restorations.


### Indication

- Used for permanent cementation of veneer made from materials of porcelain/ceramic and composites.

### Advantages

- 5 shades and 3 opacity options, satisfy all needs of esthetics, high translucency and shades covering.
- Light curing ensures colour stability.
- Impressive radiopacity, X-ray radiation resistance could be over 400%.
- High bonding strength, low water absorption and solubility, and impressive marginal fitness.
- Additional try-in gel ensures appropriate cement shade to fit in natural teeth.

### Shade



### Packaging

2.5g\*1 Syringe  
10 Dispensing Tips

## TopCEM Try-In Gel

### Veneer Try-in Gel

Delivers accurate shade matching to TopCEM veneer cement.

### Indication

- Used to evaluate the shade and the matching of the restoration prior to permanent cementation.

### Advantages

- Try-in pastes offer excellent color match to final cured cement.

### Packaging

2g\*1 Syringe  
10 Dispensing Tips



400%  
High Radiopacity



# TopCEM Veneer System

## Veneer Cementation Solution

### Indication

TopCEM Veneer System is designed for permanent cementation of veneer made from materials of porcelain/ceramic and composites.


### Packaging

6ml*1 bottle	HugeBond Universal FliPro Light Cure Dental Adhesive
2.5g*5 syringes	TopCEM-Veneer Light Cure Veneer Cement
1.5g*5 syringes	TopCEM-Try In Veneer Try-In Gel
5ml*1 bottle	TopCEM Ceramic Primer Ceramic Coupling Agent
3ml*1 syringe	P-Etchant Phosphoric Acid Etching Gel
1.2ml*1 syringe	HF-Etchant Hydrofluoric Acid Etching Gel
Accessories	Dispensing tips, mixing pad, micro-brushes

### 01 Veneer Try In



Apply Try-in gel to veneer inner face.



Carefully seat the veneer without pressure.




Rinse the veneer thoroughly, then dry with air.

### 02 Pretreatment of Veneer Restoration


• Porcelain/ glass ceramic veneers



Apply a layer of HF-Etchant, wait for 60sec.



Rinse thoroughly for 10sec, then dry with air.




Apply two layers of TopCEM Ceramic Primer to veneer inner face, and wait for 60sec.



Dry the restoration with air.

• Composite resin veneers




Apply HugeBond Universal FliPro and wait for 20sec.



Dry the restoration with air for 10sec.  
\*DO NOT LIGHT CURE.

### 03 Pretreatment of Teeth




Rinse the teeth thoroughly, then dry with air.




Isolate the teeth with rubber dam.




Apply P-Etchant to teeth surface.



Rinse the etched teeth for 10sec, then dry with air.




Apply HugeBond Universal FliPro to the teeth evenly, and rub it in for 20sec to ensure a shiny surface,




Dry with air for 10sec to remove the solvent.  
\*DO NOT LIGHT CURE.


### 04 Veneer Cementation



Apply TopCEM-Veneer Cement.




Seat the veneer in position.



Tack cure excess cement in the margin for 2sec.




Remove and clean the excess cement.



Apply glycerin gel to the margins before final curing.



Final cure the cement for 20sec.



Clean up and polish.



Apply Pro Shield Fluor Protector to the teeth for postoperative protection.





# TopCEM-RMGI

## Resin Modified Glass Ionomer Cement-Luting

Meet all your expectations for a straightforward, flexible and predictable cementation.

### Indication

- Cementation of porcelain-fused-to-metal (PFM) crowns and bridges
- Cementation of metal crowns, inlays and onlays
- All-zirconia strengthened core ceramic restorations
- Cementation of resin inlays, onlays, crowns and bridges
- Final cementation PFMS, metal crowns or all-zirconia strengthened core ceramic restorations to implant abutments

#### Simplified Operation

- No more need for adhesives, primers or conditioners.
- Tack curing ensures a quick and easy excess cement clean-up.
- Ideal consistency created for accurate placement.

#### Functional Result

- Continuous fluoride release prevents secondary caries.

#### Reliable Cementation

- Low film thickness for a marginal fit.
- High strengths delivers an ideal long-lasting cementation performance.

### Packaging

8g\*1 Syringe

### Technical Parameter

Mixing Ratio (base:catalyst)	1:1
Working Time (minimum)	1min 30sec
Intra-oral Setting Time (maximum, from initial mix)	4min 30sec

Over **20MPa**  
Reliable Bond to Zirconia



# TopCEM Vigor SA

## Self Adhesive Resin Cement

An innovation simplifies your path to restoration success.


### Indication

- Cementation of all-ceramic, composite, and metal crowns, bridges, inlays, onlays, dental posts; 2-3-unit Maryland bridges and 3-unit inlay/onlay bridges (contraindicated for patients with bruxism or periodontitis).
- Cementation of fiber posts.
- Cementation of all-ceramic, composite, and metal restorations on implant abutments.

### Advantages

- Provides a time-saving one-step cementation solution, no need of adhesive
- Ultra-low film thickness about 8µm, ensures strong luting result
- High flexural strength over 100Mpa
- Low water absorption and solubility, delivers a long-lasting restorative result.
- Perfect viscosity for accurate placement, without runny mess in patients mouth
- 7 shades and 3 opacities provide aesthetics tooth-like restorative result.


### Shade



### Packaging

8g\* 1 Syringe

### Bonding Strength(Light Curing)/ MPa



### Technical Parameter

Flexural Strength/ MPa	113.6 MPa
Film Thickness/ µm	8.4 µm
Water Absorption/ µg/mm³	27.8 µg/mm³
Solubility/ µg/mm³	0.6 µg/mm³
Working Time	≥2min
Curing Time (self curing)	2min30sec-3min10sec





TrusFIL

Universal Composite Restorative

Stands out among restoratives due to its unexpected strength and polish retention.

Indication

- Direct anterior and posterior restorations (including occlusal surfaces);
- Core Build-ups;
- Splinting;
- Indirect restorations including inlays, onlays and veneers.

Impressive Longevity

- High flexural strength, low water absorption and solubility.
- Low polymerization shrinkage effectively prevents secondary caries.
- Zirconia filler enhances compressive strength.


Easy Operation

- Shape easily, not sticky.
- Over 400% of high radiopacity for easier identification at recall.


Natrual Aesthetics

- 3 opacity and 10 shades provide natural tooth-like effect.
- High polishability and long-run gloss.

Shade



* Water absorption/ $\mu\text{g}/\text{mm}^3$	17/ $\mu\text{g}/\text{mm}^3$
* Water solubility/ $\mu\text{g}/\text{mm}^3$	0.7/ $\mu\text{g}/\text{mm}^3$
* Polymerization shrinkage/%	1.06%
* Depth of curing (A3 shade) /mm	2.49mm




Packaging

4g\*1 syringe

\* High compressive strength (321 MPa) makes composite strong enough to be placed in any class of cavity.

\* Endures load in stress bearing areas due to high flexural strength (140 MPa).

Class II restoration



\* Photo courtesy of Dr. Cindy.

Class IV restoration





# TrusFIL-Flow

Flowable Composite Restorative

Features multi-shade aesthetics, satisfying viscosity, and self-leveling technology.

Indication


- Class III and V restorations.
- Restoration of minimally invasive cavity preparations (including small, non stress-bearing occlusal restorations).
- Base/liner under direct restorations.
- Repair of small defects in indirect restorations.
- Pit and fissure sealant.
- Undercut blockout.
- Repair of resin and acrylic temporary restorations.

Advantages

- Impressive physical and chemical properties guarantee durability.
- Superior flexural strength .
- Good rheological makes it easy to control accurately and be self-leveling.
- Multiple shades provide an ideal aesthetic effect.
- Low polymerization shrinkage.
- Over 300% X-ray radiopacity .


Flexural strength/Mpa	130Mpa
Water absorption/μg/mm³	19/μg/mm³
Water solubility/μg/mm³	0.48/μg/mm³
Depth of curing (A3 shade) /mm	2.9mm

Shade



Specification

2g\*1 syringe+10 dispensing tips



130MPa  
Impressive Flexural Strength



# TrusFIL X-Blend

Universal Composite Restorative

ONE syringe for ALL ranges of teeth shades

## Indication


- Direct anterior and posterior restorations (including occlusal surfaces);
- Core Build-ups;
- Splinting;
- Indirect restorations including inlays, onlays and veneers.

## Advantages

- No wasted shades with only one TrusFIL-Universal shade, changing your restoration routine.
- Take the guesswork out of finding the accurate shades matching patients’ natural teeth.
- Exceptional longevity with low shrinkage of only 1.06%.
- Impressive physical properties due to TrusFIL Nanofiller technology.
- 3.5 mm depth of cure – saves time.
- Easy to polish, highly stain resistant with a high gloss.


## Specification

4g\*1 Syringe



## TrusFIL X-Blend Technology

- Utilizing special fillers and innovative TrusFIL Technology, TrusFIL X-blend exhibits succeeds in covering all shades in only ONE universal shade composite with the Chameleon Effect.



## Properties

Flexural Strength/Mpa	138 Mpa
Polymerization shrinkage/ %	1.06%
Elastic Modulus/ GPa	11.5 GPa
Water Absorption/ µg/mm³	18.2 µg/mm³
Solubility/ µg/mm³	0.5 µg/mm³
Depth of Curing/ mm	3.5 mm

\* All data from HUGE lab





# TrusFIL Restoration System

## Universal Composite Restoration Solution

### Indication

TrusFIL Restorative System is a highly aligned direct restorative system indicated for:

- Direct anterior and posterior restorations (including occlusal surfaces);
- Core Build-ups;
- Splinting;
- Indirect restorations including inlays, onlays and veneers.

### Packaging

6ml*1 bottle	HugeBond Universal FliPro Light Cure Dental Adhesive
4g*4 syringes	TrusFIL Universal Composite Restorative
2g*2 syringes	TrusFIL Flow Flowable Composite Restorative
3ml*1 syringe	P-Etchant Phosphoric Acid Etching Gel
Accessories	Dispensing tips, mixing pad, micro-brushes





## UltraCore

### Core Build-up Material

A real-dentin-like core build-up material, brings you utmost restorative results.


### Indication

- Core build-up.
- Root post cementation.
- Permanent crown&bridge cementation.

### Advantages


- Superior compressive and flexural strength.
- Cuts like natural dentin and matches natural tooth structure.
- Excellent viscosity guarantees optimal adaptation to the cavity walls and root posts.
- Ideal Ultra-thin film thickness (12.7μm) for cementing root posts.
- 7-shade options ensure natural tooth look.
- Dual curing techniques provide easier application in the root canal.

### Packaging




### Optimal Material for Core Build-ups


Compressive Strength/MPa



Flexural Modulus/ GPa



### Root post cementation



### Core Build-up





## Proseal F

### Enamel Coating Resin

A safe and painless way of protecting teeth from tooth decay.


### Indications

- Deep pits and fissures on the occlusal surface, buccal surface and lingual or palatal surface, especially where a probe can be inserted or stuck (including suspected caries)
- The contralateral tooth has caries or tends to caries
- Tooth eruption reaches the occlusal surface or the pits and fissures on the crown surface are completely exposed in the oral cavity


### Advantages

- Upgraded tip design enables accurate placement, minimizes leakage, and virtually eliminates bubbles.
- Low viscosity allows easy flow into pits and fissures
- Releases fluoride stably to prevent cavities
- Changes from pink to white after light curing
- Offers excellent stain resistance
- No pungent odor

### Parameters



\* Data from HUGE laboratory



\* Data from HUGE laboratory

### Operation Guide


as Pit & Fissure Sealant



Prepare the tooth (Rinse and dry).



Etch the pit and fissures with P-Etchant Phosphoric Acid Etching Gel.



Rinse and dry.



Apply Proseal F to the tooth.



Light cure until Proseal F turns to opaque white colour.



Finish the prevention.

### Packaging

1.5g\*1 Syringe  
10 Dispensing Tips





## Pro Shield

### Fluor Protector

An innovative varnish system with a stable suspension of sodium fluoride.

### Indication

- Treating hypersensitive teeth



### Advantages

- High-efficiency fluoride release
- Safe natural resin formula
- Stable and even fluoride suspension
- Effective relief of tooth hypersensitivity
- Pleasant flavor for better patient experience
- Clinician-friendly packaging design
- Light yellow color from natural resin provides visual control

### Packaging

- 10g/tube
- 0.4g\*50/box

### High Fluoride Release Capability Suitable for Kids and Adults



\*Data from HUGE laboratory


### Excellent Operation Experience with Pro shield



1. Make sure the areas to be treated are clean




2. Remove excess saliva with triple syringe or cotton pellets




3. For 10g tube:



4. For 0.4g unit-dose:



5. Use the microbrush to apply a thin layer to the tooth surface



6. Forms a thin film on the tooth surface

**22,600 ppm**  
**Sodium Fluoride**



# Pediatric Dental Solution

## Caring for Children's Dental Health To benefit you and your young patients.

HUGE has a professional pediatric dental system. The range includes preventive, minimally invasive treatments and restorative products. All products have shown effectiveness and safety in clinical applications, helping children maintain good dental conditions and effectively treating children's dental diseases.

### Why choose glass ionomer cement materials for children?

- Continuous and long-lasting fluoride release to prevent children's dental cavities
- Easy to operate and simple to apply clinically
- Excellent bonding strength to metal restorations. Can be used to bond prefabricated crowns

### Prevention

To keep teeth healthy, small amounts of fluoride are needed during all stages of tooth growth. Especially for children, the prevention of dental cavities is very important. Several dental organizations and research reports recommend the professional use of fluoride restorative materials.

### Diagnosis

**Proseal F**  
Enamel Coating Resin



**Pro Shield**  
Fluor Protector



### Restoration

**TrusFIL**  
Unviersal Composite Restorative



**TrusFIL-Flow**  
Flowable Composite Restorative



### Protection

**Glass Ionomer Cement**  
Universal Enhanced Glass Ionomer System  
Luting I · Filling I · LuFill HS



**HugeBond Universal FliPro**  
Light Cure Dental Adhesive









# PERFIT Putty

## Elastomeric Impression Material

### Indications

- For crown, bridge, inlay, and onlay impressions.
- Recommended for use with wash material in two-step impression technique.

### Advantages

- Good linear dimensional stability.
- Excellent recovery after deformation.
- High tear strength and tensile strength.

### Technical Parameters


Technical Parameter	Putty	
	Normal Set	Fast Set
Type	Normal Set	Fast Set
Mixing Time	30s	30s
Working Time	1min30s	1min15s
Setting Time in Oral Cavity	3min00s	2min30s
Linear Dimensional Change (24 h)	>99.5%	
Recovery From Deformation	<0.2%	
Shore A Hardness (24 h)	≈67	
Color	Purple	
Shelf Life	3 Years	

\*Data from HUGE laboratory



**Packaging**

400g*1	Base
400g*1	Catalyst
2	Measuring Spoons



**Packaging**

400g*1	Putty Base
400g*1	Putty Catalyst
2	Measuring Spoons
50ml*1	Light Body Cartridge



# PERFIT Heavy Body



## Elastomeric Impression Material

### Indications

- For crown, bridge, inlay, onlay, implant, orthodontic, and full denture impressions.
- Recommended for use with wash material in one-step or two-step impression technique.

### Advantages

- Easy and efficient operation.
- Good mixing quality, avoid bubbles.
- High precision and dimensional stability.
- High tear strength and tensile strength.
- Especially suitable for border molding.



Heavy Body/Dispenser-Mix



Heavy Body/Auto-Mix

### Impression Technique

1. Double Step Technique  
Recommended for all clinical restorations, crowns, bridges, inlays, onlays and implant impressions
2. Single Step technique  
Recommended for implant impression



Packaging  
380ml\*1 Cartridge



Packaging  
50ml\*2 Cartridges



Packaging  
50ml\*2 Cartridges  
6 Mixers

### Technical Parameters

Technical Parameter	Heavy Body/Dispenser-Mix		Heavy Body/Auto-Mix
Type	Normal Set	Normal Set	Fast Set
Working Time	2min00s	1min30s	1min30s
Setting Time in Oral Cavity	3min30s	3min00s	3min00s
Linear Dimensional Change (24 h)	>99.5%		
Recovery From Deformation	<0.2%		
Shore A Hardness (24 h)	≈68		
Color	Blue		Green
Shelf Life	3 Years		

\*Data from HUGE laboratory



# PERFIT Regular Body-Transparent

## Elastomeric Impression Material

### Indications

- For Clinical Application:
  1. Two-step impression technique with wash materials
  2. Light-curing guide for
    - a. Mock up restoration
    - b. Bonding of orthodontic brackets
    - c. Occlusal stamp for composite restorations
- For Lab Application:
  - Composite injection moulding technique



### Advantages



Extreme Transparency



No bubbles



Ideal Rigidity



High Aesthetic Results

### Technical Parameters

Technical Parameter	Regular Body
Working Time	1min30s
Setting Time in Oral Cavity	3min00s
Shore A Hardness (24 h)	62
Color	Transparent
Shelf Life	3 Years


\*Data from HUGE laboratory



**Packaging**  
50ml\*2 Cartridges

### Comparison


PERFIT Regular Body in comparison to competitors: From the comparison, it can be clearly seen that HUGE Regular Body has better transparency. And it rarely produces air bubbles and has high aesthetic effect.




### Clinical Application Case

The matrix created with PERFIT Regular Body achieves excellent aesthetics. In a short period of time, it can simplify clinical procedures and increase operational efficiency.


### Flowable Injection Technology for Aesthetic Restoration




1. Prepare teeth and make the plaster model




2. Finish restoration with wax on plaster model




3. Apply Regular Body to the plaster model



4. Use Regular Body to make transparent matrix



5. Fill the space with TrusFIL-Flow (Flowable Composite Restorative)



6. Light cure the resin materials



7. Finish restoration after polishing



# PERFIT Light Body



## Elastomeric Impression Material

### Indications

- For crown, bridge, inlay, onlay, implant, orthodontic, and full denture impressions.
- Recommended for use with tray material in two-step or one-step (simultaneous) impression techniques.

### Advantages

- Extraordinary flow property and stability.
- High hydrophilicity and outstanding thixotropy.
- Strong ability to reproduce details.



### Technical Parameters

Technical Parameter	Light Body	
Type	Normal Set	Fast Set
Working Time	2min15s	1min30s
Setting Time In Oral Cavity	3min30s	2min30s
Linear Dimensional Change (24 h)	>99.5%	
Recovery From Deformation	<0.2%	
Shore A Hardness (24 h)	≈56	
Color	Orange	
Shelf Life	3 Years	

\*Data from HUGE laboratory



**Packaging**  
50ml\*2 Cartridges



**Packaging**  
50ml\*2 Cartridges  
12 Mixers  
12 Oral Tips


### Clinical Application



01. Choose PERFIT Light Body wash material



02. Apply the material on the initial impression and in the patient's mouth



03. Complete accurate impression



Impression results  
can be stored for  
**28 days**



# PERFIT Super Light Body

## Elastomeric Impression Material

### Indications

- Especially recommended for full denture impressions.
- Recommended for use with tray material in two-step impression technique.

### Advantages

- Higher flowability than Light Body with excellent stability.
- High tear strength and tensile strength.
- Superior accuracy for precise detail reproduction.



### Technical Parameters

Technical Parameter	Super Light Body
Type	Normal Set
Working Time	2min15s
Setting Time In Oral Cavity	3min30s
Linear Dimensional Change (24 h)	>99.5%
Recovery From Deformation	<0.2%
Shore A Hardness (24 h)	≈56
Color	Dark Green
Shelf Life	3 Years

\*Data from HUGE laboratory




**Packaging**  
50ml\*2 Cartridges



**Packaging**  
50ml\*2 Cartridges  
12 Mixers  
12 Oral Tips

### Clinical Application





# PERFIT Heavy Body & Light Body Fast Set



## Specific Applications Solutions:

Compared with traditional two-step impression technique, Heavy Body&Light body benefits from its simple operation and excellent dimensional stability. We recommend it to be used as an impression tool for implant restoration and full denture restoration.



## Unique Benefits:


- High hydrophilicity for precise and reliable impressions.
- High tear strength and tensile strength.
- Excellent linear dimensional stability help to maintain impression accuracy for up to 4 weeks.
- Good bio-compatibilitysafe to use for patients




## Dental Advisor Evaluators' Comments

"The ease of loading and consistency of this impression material are very good."  
"Great detail and handling. Excellent hydrophilic property."  
"The color of the light body made it easy to read margin details."  
"I loved the set-up time!"  
"Pretty good setting time"  
"Very clean impressions, easy to express, hydrophilic and worked with our impression tips that seemed to give less bubbles"

"GREAT DETAIL AND HANDLING. EXCELLENT HYDROPHILIC PROPERTY."



Compared to Competitive Products:





# Impression Kit

## Advantages



Total Solution  
Provider



Rich Customized  
Options



Reliable  
Quality

## Versatile Solutions

Widely used in on-line oral treatment like invisible orthodontic, teeth beauty and sports mouthguard.

## HUGE Standard Package

- Putty Base 25g\*4+ Putty Catalyst 25g\*4
- Upper trays\*2 + Lower trays\*2
- Instruction for Use\*1



## Customized Options

- 01** Full Solution Provider  
Integrated product supply
- 02** Rich Customized Options  
Diversified product configuration combination
- 03** Reliable Quality  
Same material with HUGE Perfit





# PERFIT Bite Registration

## A-Silicone for Bite Registration


90 Shore A  
High Hardness

### Indications

- For making bite registration.

### Advantages

- Accurate bite record with high hardness(90 A).
- High rigidity and easy trimming.
- Excellent dimensional stability.
- Whipped-cream-like consistency.
- Scannable for digital workflows.



### Technical Parameters

Technical Parameter	Bite Registration	
Type	Normal Set	Fast Set
Working Time	30s	12s
Setting Time In Oral Cavity	1min30s	45s
Shore A Hardness (24 h)	90	
Color	Light Blue	
Shelf Life	2 Years	

\*Data from HUGE laboratory



Packaging  
50ml\*2 Cartridges



Packaging  
50ml\*2 Cartridges  
6 Mixers



# PERFIT Tray Adhesive

## Impression Accessory

### Indication

-To offer reliable adhesion between impression tray and silicone impression material.

### Advantages

- Suitable for A-Silicone and C-Silicone
- Easy to apply and handle
- Short setting time after application
- Blue colour visual control

### Drying Time

3 minutes



### Packaging

10ml/bottle

### Operation Guide




01




02

Apply a uniform layer of PERFIT Tray Adhesive onto the tray, dry for at least 3 minutes, then apply the silicone impression material.








# Technical Test Results

Precise · Professional · Efficient


## Waterdrop Test ▼



Contact angles of Light body with different brands

Excellent hydrophilicity helps obtain the details quickly and reduce the risk of bubble formation,thus guarantee a high level of precision.

Hydrophilicity test in 2 secs



\*According to HUGE laboratory test results



## Drip Test ▼

Outstanding thixotropy ensures a better clinical experience.




High flowability allows to reach the less accessible spaces of the gingival sulcus

## Flowability Test ▼

















\*According to HUGE laboratory test results



\*According to HUGE laboratory test results



# Silicone Impression Material

	Dispensing system	Standard Package	Accessory	
Putty		400g*2		
Light Body		50ml*2		
Super Light Body		50ml*2		
Heavy Body Dispenser-Mix		50ml*2		
Bite Registration		50ml*2		
Regular Body		50ml*2		
Heavy Body Auto-Mix		380ml*1		

## Stable

No deformation up to 4 weeks

Due to excellent linear dimensional stability, impressions made from PERFIT are very stable and can be stored for up to 4 weeks without deformation.  
Excellent silicone impression material will keep the impression dimensional change rate < 0.2%. This is an important condition for obtaining accurate impressions. As a result, the technician can create a more precise restoration based on the impression.







# Implant Impression Solution

HUGE is dedicated to implant restoration and has a complete product line from clinical products to laboratory materials. Precise impression-taking is one of the important steps in the process of implant restoration. A tight fit between the impression material and the impression post is required to achieve an accurate recording of the seating orientation of a restoration, whether an open tray impression technique or a closed tray impression technique.



## PERFIT IMPLANT IMPRESSION WORKFLOW

### Precise impression-taking



**STEP 1**  
Apply the **PERFIT Putty/Heavy Body**



**STEP 2**  
Apply the **PERFIT Light Body**




**STEP 3**  
Take impression with **PERFIT Impression materials**



**STEP 6**  
Finish the bite registration



**STEP 5**  
Apply the **PERFIT A Silicone for Bite Registration**



**STEP 4**  
Finish the implant impression

### Temporary restoration



**STEP 1**  
Make transparent matrix with-  
**PERFIT Regular Body/ PERFIT Bite Registration-Transparent**



**STEP 2**  
Apply **temporary crown and bridge resin material** in matrix



**STEP 3**  
Apply the matrix and material in the patient's mouth



# Full Denture Impression Solution

Full denture restoration requires high accuracy of impressions. Elastic impression materials are suitable for edge shaping, which is beneficial to better capture the details of the patient's mouth. Applying bite registration materials to record the horizontal bite relationship is also a necessary process for full denture restoration. HUGE provides dentists with excellent products and services to complete full denture restoration.



## PERFIT FULL DENTURE IMPRESSION WORKFLOW

### Precise impression-taking



### Full denture restoration

#### HUGE Soft Denture Relining Materials Kit

For permanently soft relining of partial or full dentures

Soft Relining material: 50ml Cartridge

Primer: 10ml /Bottle

Varnish: 5ml base, 5ml catalyst /Bottle



COMING SOON





# GLASS IONOMER CEMENT







# Glass Ionomer Cement

- Excellent biocompatibility and good tissue response
- Easier to clean than resin-bonded materials
- No etching or bonding required for fast, easy and simple application
- Ideal products for pediatric, geriatric and special needs patients

Indications		Luting I	Filling I	LuFill HS
Inlay/Onlay/ Crown/Bridge	Ceramic/Glass Ceramic	●		●
	Oxide Ceramic	●		●
	Metal/Metal-Based	●		●
	Composite Resin			
Post/Screw	Metal	●		●
	Oxide Ceramic	●		●
Orthodontic Bands		●		●
Primary Teeth			●	●
Class III, V and Limited Class I Cavities			●	●
Base Lining			●	●
Core Build-Up			●	●
Atraumatic Restorative Treatment			●	●
Cavity of Unstressed Area			●	●



# Glass Ionomer Cement

## Luting I

A truly biocompatible, quick-set, contemporary luting glass ionomer cement that ensures reliable and long-lasting restorations.

### Advantages

- Chemically bond to tooth structure for excellent edge sealing
- Excellent biocompatibility and good tissue response
- Optimum thixotropy and ideal film thickness (12µm)
- Excellent bonding strength (>6MPa)
- Easier to clean than resin-bonded materials

### Technical Parameters

Luting I	
Mixing Time	45s
Working time	2min10s
Net setting time	3min05s
Film thickness	12µm
Compressive strength	>150 MPa

\*Data from HUGE laboratory

### Packaging

- 30g powder + 25g liquid +1 measuring spoon + 1 mixing pad
- 10g powder + 10g liquid +1 measuring spoon + 1 mixing pad

### Product color after curing



14±2 µm  
Ideal Particle Size



# Glass Ionomer Cement

## Filling I

A true self-adhesive, bulk filling glass ionomer cement that ensures reliable and long-lasting restorations.

## Advantages

- Easy to mix and provides good operation experience
- High compressive strength (>180MPa)
- Ideal particle size distribution
- Continuous fluoride releasing
- Unique color system provides good aesthetic effect
- X-ray radiopacity

## Technical Parameters

Filling I	
Mixing Time	45s
Working time	2min20s
Net setting time	2min45s
Compressive strength	>180 MPa

\*Data from HUGE laboratory

## Packaging

- 15g powder + 12g liquid +1 measuring spoon + 1 mixing pad
- 5g powder + 5g liquid +1 measuring spoon + 1 mixing pad

Available in Shade A1,A2,A3

\*HUGE unique color system



7.0  $\mu\text{m}$   
Average  
Particle Size



# Glass Ionomer Cement

## LuFill HS

A truly self-adhesive, **high-strength** universal glass ionomer cement that ensures a reliable, long-lasting restoration.

### Advantages

- Easy to mix and provides good operation experience
- High compressive strength (>200MPa)
- Low stickiness to appliance
- Continuous fluoride releasing
- Universal Vita color system provides good aesthetic effect
- X-ray radiopacity

### Technical Parameters

Technical Date	LuFill HS	
	For Luting	For Filling
Mixing Time	45s	45s
Working time	2min10s	3min00s
Net setting time	3min15s	2min45s
Film thickness	13µm	/
Compressive strength	/	>200MPa

\*Data from HUGE laboratory

### Packaging


- 30g powder + 25g liquid +1 measuring spoon + 1 mixing pad
- 10g powder + 10g liquid +1 measuring spoon + 1 mixing pad

Available in Shade A2,A3,A3.5

\*Vita color system



### Technical Parameters



>200MPa  
Compressive Strength



# GLASS IONOMER CEMENT CAPSULE

## Luting I • Filling I


HUGE Glass Ionomer Cement Capsules, the new packaging is easier to mix and convenient to operate. It still maintains the excellent performance of traditional GIC. It is available in Filling and Luting.

- Universal Application**  
It is available in Luting I and Filling I and has clinically proven performance.
- Fast and simple**  
Capsule packaging is easy and quick to use.
- Ideal consistency**  
Perfect mix ratios optimize flow, cavity fit, reduce voids and ensure strongest chemical adhesion.
- High compressive strength(>180MPa)**  
The precise ratio of powder and liquid ensures maximum physical properties, resulting in a stronger, more durable restoration.
- Convenient and consistent**  
The same predictable result each and every time.

## Packaging

Luting I Capsules : 50 capsules/Box  
Filling I Capsules : 50 capsules/Box


## Operation guide



1. Shake or tap to loosen powder.



2. Depress plunger, hold down 2 seconds.




3. Mix for 10 seconds.



4. Two clicks to prime capsule then syringe slowly.

COMING SOON





# INTRAOURAL SCANNER





# Realistic color scanning



to help show a more detailed 3D image of the teeth

A powder-free intraoral scanning device can display the teeth details in a high-definition natural view based on video dynamic intraoral scanning.

Starting 3D light structure scanning technology, high-speed and anti-shaking performance

AI technology makes you enjoy the smart experience and smooth scanning.

Intelligent tracking technology ensures fast recovery scanning of any oral position.




# Multiple scanning heads, Infinite options

Multi functionality, safe and reliable, universal data format


## Restoration

Accurate and clear shoulder edge guarantees more efficient design.  
True color scanning can effectively distinct the gums and tooth structure.  
Fast scanning results permit immediate restoration.




## Orthodontics

The full dentition scan has high accuracy and restores the true shape of the dental arch.  
Immediate treatment plan, low time-cost.



## Implant

Fast scanning speed, accurate scanning metal positioning rod.  
Great scanning depth, available 3mm data under cuff.  
Improved experience, no repeated scanning.






# Software Highlights


# Transmission in One-Click

Coming Soon




## Undercut Area Analysis

It allows to check the undercut at different perspectives. The undercut area will be marked with gradient color.




## Align

Recalculate the current occlusion relationship.



## Distance


It allows to calculate the occlusal distance and show the result, and to separate the upper and lower jaw to check the distance.



## Lock Area

It is used to protect the area that you do not want to change during the scanning process.

Transmit your scanning data within the software, it's high-efficient and time-saving. The real-time data process technology strengthens treatment performance.



ANALYTICS



REAL-TIME



USER-FRIENDLY  
HMI



# Accurate, Precise, and Fast

Bring an unprecedented digital scanning experience



AI technology improves intelligent scanning experience.



Clear cervical margin, vivid color and high precision.



A variety of scanning tips to meet various complex cases.



Fast scanning speed, easy to use.




Portable and micro device, easy to operate.



Powder free for true color scanning.

HUGE Scanner	Traditional impression
<p>✓ <b>Short treatment period</b></p> <p>Data can be obtained directly, 1 minute to take full mouth model.</p>	<p><b>Long treatment cycle</b></p> <p>It takes 5-10 minutes for impression and usually takes 1-2 weeks to complete the dental restoration. The rework rate is high.</p>
<p>✓ <b>High data accuracy</b></p> <p>High accurate scanning data, and it is convenient for date storage and trace.</p>	<p><b>Low data accuracy</b></p> <p>Indirectly data access through multiple processes, e.g. impression, reversal, and scanning, which reduces the accuracy of data.</p>
<p>✓ <b>Comfortable</b></p> <p>Short mouth opening time, no abnormal sensation in the mouth, and available to be paused at any time.</p>	<p><b>uncomfortable</b></p> <p>Long patient's mouth opening time impacts impression taking experience, and the abnormal sensation will lead to nausea and vomiting.</p>
<p>✓ <b>Informative</b></p> <p>The visual digital model can help analyze and simulate intraoral changes in time, such as tooth wear, gingival recession, etc.</p>	<p><b>less information</b></p> <p>It is impossible to observe the dynamic changes based on the single static model data.</p>


HUGE Smart	Specification
Technology	Staring Scanning with Coded Structure Light
Anti Fog System	Intelligent Heating
Accuracy	≤ 7μm
Precision	≤ 5μm
Handpiece Weight	138g (192g with cable)
Handpiece Dimension	221×27×25mm
Handpiece Enclosure	High-Strength Aviation Aluminum Alloy
Scan Field	18×16mm
Scan Depth	0-20mm
True Color	Full HD
Tips	6
Tip Autoclave	50-100 times
Scanning Control	Button
Motion Sensing	Support
Cooling	Turbo Fan
Calibrator	Automatic
Data Interface	USB3.0
Cable Length	1.85m





# HUGE provides scenario-based solutions

HUGE is dedicated to clinical aesthetic restoration, providing dentists with a comprehensive portfolio of high-quality materials and professional clinical solutions.



**Implant Restoration Solutions**

HUGE focuses on implant solutions. From initial impressions to temporary restorations and final restorations, HUGE will be your most reliable partner.



**Prevention and caries infiltration**

Effective prevention is a priority, especially in the field of dentistry. HUGE has been focusing on the prevention of caries, especially the early protection of children's teeth.

**We are here for you. HUGE**

- All products of HUGE brand are carefully produced, strictly controlled and certified.
- Please believe that we have the ability to provide you with safe, reliable and high-quality clinical products.

Welcome to visit HUGE official website and store:

[www.hugedental.com](http://www.hugedental.com)


FDA Approved CE 1639 ISO 13485:2016



# HUGE promises professional training

As a partner of dentists, HUGE aims to support you in all aspects, including product solutions and practice-related professional training activities.

With a professional R&D and clinical team, HUGE has carried out a variety of themed activities for clinical applications.



## All Knowledge. All Free.

HUGE provides a variety of educational videos and live online courses, introducing HUGE's products and scientific background, and providing professional operation guidance.


Contact us for more details if you want to know.



HUGE Regular Body(Transparent) Operation Guide 2



HUGE Regular Body(Transparent) Operation Guide 1



Huge Clear Aligner OEM Service



HUGE Bite Registration Silicone Operation Guide