

## CASTING WAX Instructions for use

Version: 01  
Issuance Date: Dec. 08, 2021

【Product name】 Casting Wax

【Product composition】 Wax

【Application】 To fabricate male mold for oral hard and soft tissue or restoration model

【Class】 Hard

【Product performance】

Items	Technical parameters
Melting point	108.4℃
Ignition residue at 500℃	0.051%

【Operation】










According to the routine operation manner of the dental technician

【Storage】 Indoor storage at room temperature

【Warning】 None

【Precaution】 Notification of any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority in your country.

【Symbol】

 --- Batch code     --- Manufacturer ---  --- Country of manufacture  
 --- Medical device     --- Caution -----  --- Unique device identifier  
 --- Consult instructions for use  
 --- Date of manufacture  
 --- Authorized representative in the European Community

**Rx** ONLY .... Caution: U.S.Federal Law restricts this device to sale by or on the order of a dentist or professionals.

 Shandong Huge Dental Material Corporation  
No. 68 Shanhai Road, Donggang District, Rizhao City, Shandong  
Province, 276800, P.R. China  
Tel: +86 633 2277 268  
marketing@hugedental.com www.hugedental.com

 MedNet EC-REP GmbH  
Borkstrasse 10, 48163 Muenster, Germany

